## Revolix HTL


Hybrid Thulium Laser Family



Versatile performance, no compromises

**Revolix HTL** is the powerful All-in-One Hybrid Thulium Laser based on proven **Revolix** Technology, customized for superior stone and tissue surgery.





## **Advantages**



### Stones

Superior stone dusting and fragmentation compared to Holmium systems at equivalent settings.

Simplified treatment and increased stone control due to negligible retropulsion effect.

Faster and more precise disintegration of concrement due to high frequencies.



#### Tissue

Superior pulsed tissue enucleation performance delivers *RevoLix* proven haemostasis with significantly lower carbonization effect.

High power pulsed mode leads to equivalent results compared to HoLEP and ensures effective vaporization of tissue.



### **Economic**

Revolutionary Hybrid Thulium Technology supports all laser techniques and allows faster procedures.

RevoLix HTL devices are compatible with both LISA disposable and reusable fibers.

Green Engineering: Lower energy consumption compared to Holmium Laser devices.

# Strengthen your performance



	RevoLix HTL	RevoLix HTL eco	RevoLix HTL cw		
Wavelength	2,013 nm	2,013 nm	2,013 nm		
CW mode					
Emission mode	continuous wave	continuous wave	continuous wave		
Power	5 - 150 W, adjustable in 1 W steps	5 - 75 W, adjustable in 1 W steps	5 - 150 W, adjustable in 1 W steps		
PULSE mode (QCW)					
Emission mode	pulsed	pulsed	-		
Average Power	5 - 150 W	5 - 75 W	-		
Pulse peak power	max. 1,300 W, adjustable	max. 1,300 W, adjustable	-		
Pulse repetition rate	5 - 300 Hz	5 - 300 Hz	-		
Pulse duration	200 - 4,750 μs	200 - 4,750 μs	-		
Pulse energy	max. 4.5 J, adjustable	max. 4.5 J, adjustable	-		
Aiming beam	532 nm (green), 3 mW (adjustable)	532 nm (green), 3 mW (adjustable)	532 nm (green), 3 mW (adjustable)		
Mains supply	200 V - 240 V, 50/60 Hz, 10 A (1~, N, PE) or 110 V - 115 V, 50/60 Hz, 20 A (1~, N, PE), 20 A (1, N, PE), automatic range selection		200 V - 240 V, 50/60 Hz, 10 A (1~, N, PE) or 110 V - 115 V, 50/60 Hz, 20 A (1~, N, PE), automatic range selection		
Cooling system	Air cooling	Air cooling	Air cooling		
Dimensions (HxWxL)	1,025 mm x 450 mm x 740 mm	1,025 mm x 450 mm x 740 mm	1,025 mm x 450 mm x 740 mm		
Weight	108 kg	108 kg	108 kg		
Environmental conditions	15 - 28 °C / 10 - 90 % humidity	15 - 28 °C / 10 - 90 % humidity	15 - 28 °C / 10 - 90 % humidity		

## Simple and Safe

- ✓ Hybrid operation modes
- Memory function
- ✓ Harmonized Footswitch Control
- Standard Power Plug

Integrated Whisper Technology<sup>TM</sup> makes RevoLix HTL one of the quietest high power lasers in the OR.



## **LISA Laser Fibres**

Quality, Power and Precision

## **Single Use Sterile**

	<b>LithoFib-SU</b> GlideTip™ bare fibre	<b>SureFib-SU</b> GlideTip™ bare fibre	<b>FlexiFib-SU</b> GlideTip™ bare fibre	PercuFib-SU bare fibre	<b>RigiFib-SU</b> bare fibre	<b>RigiFib-SU 800</b> bare fibre	<b>RigiFib-SU 1000</b> bare fibre	<b>SideFib-SU</b> side fibre deflection 70°	
Order no. of fibre	101 503 576	101 503 513 101 503 569	101 503 387	101 503 384	101 503 289	101 503 385	101 503 386	101 503 138	
Optical Core / micron	200	272	272	365	550	800	940	550	
Length / m	2.5	2.5	2.5	2.5	2.5	2.5	2.5	3	

## **Reusable fibres**

	<b>LithoFib</b> bare fibre	<b>SureFib</b> bare fibre	<b>FlexiFib</b> bare fibre	<b>PercuFib</b> bare fibre	<b>RigiFib</b> bare fibre	<b>RigiFib 800</b> bare fibre	<b>RigiFib 1000</b> bare fibre	
Order no. of fibre	101 503 188	101 503 364	101 503 189 101 503 118	101 503 128	101 503 213	101 503 287	101 503 284	
Optical Core / micron	200	272	272	365	550	800	940	
Length / m	3	3	3 5	3	3	3	3	
FibreStrip-S Ø / mm	0.2	0.3	0.3	0.5	0.7	1.0	1.1	
Order no. of stripper	101 503 273	101 503 274	101 503 274	101 503 275	101 503 277	101 503 280	101 503 382	

Refer to the product labelling for indications, contraindications, warnings and instructions for use. Specifications are subject to change without notice.

#### Device is not for sale in U.S.A.

LISA Laser Products GmbH is a wholly owned subsidiary of OmniGuide Holdings. For more information call +49 5556 9938 0 or visit **www.lisalaser.de** 

+49 5556 9938-0

**⇒** +49 5556 9938-10

sales@lisalaser.de www.lisalaser.de

LISA Laser Products and Sphinx are registered trademarks of LISA Laser Products GmbH. All content is copyright of LISA Laser Products GmbH. © 12/2020 See www.omni-guide.com/patents for information about the coverage of LISA Laser Products by U.S. and/or foreign patent rights. To request copies use:  $0.36\,0.26\,0.00$ 



VISIBLE AND INVISIBLE LASER RADIATION AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION CLASS 4 LASER PRODUCT