A.4019L

A.4021L

A.4022L

Self · led · light





CE


Disposable **Self-LED-light** are made of transparent plastic material and supplied in individually sealed bags with LED light built-in.

LED lighting can be performed by pulling out the transparent strip. Guaranteed same light intensity for two hours usage. At the end of the examination, by lifting the tab placed on the base of the penlight it will open in two parts allowing the separate disposal of its parts (metal, batteries, plastic).

Compatible, through adapter, to systems of image aquisition (cameras) **ALREADY AVAILABLE IN HOSPITALS** (for example operating theatres laparoscopic column).

W.000160GB REV 02 03 2019

Self · *led* · light



TECHNICAL FEATURES