

Hall[®] MicroFree[™] Cordless Small Bone Power System

The NEW cordless MicroFree System is designed to make small bone procedures easier for surgeons, combining enhanced ergonomics for increased comfort, with freedom of motion for precision handling.



HALL 50° POWERED INSTRUMENTS



ASSISTARM[™] LIMB POSITIONER

T

2

TRINITY[™] INTEGRATED CAMERA SYSTEM





HALL

To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.

NEW

Unleash Your Precision.



An Easier Day In The OR.

When we asked surgeons how we could make their jobs easier during small bone cases, the answer was clear – Cordless power handpieces that provide freedom of motion without sacrificing power. It was clear these handpieces must be comfortable, ergonomic, well-balanced and deliver the same power they were accustomed to with corded systems.

Surgeons wanted a new system that delivered improvements without drawbacks – convenience without trade-offs. That is why we developed the MicroFreeTM Cordless Small Bone Power System.

It's time to cut the cord on your small bone power.

HALL® MICROFREE™ SMALL BONE POWER SYSTEM



MicroFree Pencil Grip



HALL[®] MICROFREE[™] SMALL BONE POWER SYSTEM

Cordless Precision.

MicroFree[™] went through extensive design testing with leading orthopedic surgeons¹ to create a cordless handpiece that makes it easier to perform with precision. Using their feedback, the center of gravity was adjusted to improve balance and control while reducing muscle fatigue.¹ The result was a cordless system that provides freedom of motion in the familiar pencil-grip design.

A Full Cordless Small Bone System.

One Complete System for Your Clinical Needs

The MicroFree[™] System can be used for most small bone procedures of the foot, ankle, hand and wrist. For added convenience, the sagittal saw head can be rotated in 45° increments to provide additional flexibility during cases. The MicroFree System also features the modular Mini-Driver[™] for pin/wire-driving and drilling.

SYSTEM INCLUDES:

- 5 HANDPIECES
- MICROFREE LITHIUM POWER UNIT
- MINI-DRIVER
- SMALL LITHIUM BATTERY
- CUSTOMIZABLE STERILIZATION TRAYS





Hall Mini-Driver





Small Lithium Battery

Powerful High-Performance Lithium-Ion Technology

A cordless system is only an asset if it is powerful enough to get the job done. Our high-performance lithium-ion technology provides the torque, power and performance needed.




Lithium Power Unit

HALL® MICROFREE[™] SMALL BONE POWER SYSTEM

Facility Benefits: One System That is Easy to Maintain and Upgrade

Autoclavable for Efficient and Cost-Effective Sterile Processing

Both the UL-approved batteries and handpieces of the MicroFree[™] System are autoclavable, eliminating the need for more expensive STERRAD battery sterilization. This efficient system allows the batteries, handpieces and attachments to be sterilized together in the same tray at the same time, simplifying the process, saving OR staff time and energy.



IPX 6 and 8 - Protection from Moisture

Just like our large bone Hall 50[™] handpiece, the MicroFree handpiece is rated IPX 6 and IPX 8 to help protect the system against water damage. The UL-approved hermetically sealed power units are fully autoclavable allowing the whole system to be sterilized together.

Safety – Designed to Reduce the Risk of Infection

By eliminating unsecured power cords from compromising the sterile field, the MicroFree System may reduce those potential risks that can contribute to surgical site infections (SSIs).



Combined with the Hall 50[™] System for Large Bone, the MicroFree[™] System can outfit your facility with a power solution that is completely autoclavable.

Backward and Forward Compatibility for Easy Upgrades

Often when other companies roll out new systems, they are intentionally incompatible with their old technology – forcing expensive upgrades when your existing equipment is still working. Like most Hall Systems, MicroFree is designed to be both backward and forward compatible – delivering greater long-term value for your investment.

High Quality, Flexible Service Options

As a FDA-regulated and ISO 13485 certified facility, CONMED service provides fast, cost-effective and high-quality service for every component of the MicroFree System. We also offer complete flexibility in our service solutions and multiple program levels that can be tailored to meet the specific needs of your organization's surgical schedule, inventory and budgetary requirements.

Experience the Freedom of MicroFree – Schedule an Evaluation Today

To take the next step toward discovering the comfort and convenience of cordless small bone power, schedule an evaluation by calling **1-800-237-0169** to locate your local CONMED Sales Representative.

For more information, visit ConMed.com/MicroFree.

HALL® MICROFREE™ SMALL BONE POWER SYSTEM



525 French Road, Utica, NY 13502

Ordering Information

Description	Catalog Number
NEW Hall [®] MicroFree [™] Handpieces	
Hall MicroFree High Speed Drill	PR08000SB
Hall MicroFree Medium Speed Drill	PR08100SB
Hall MicroFree Sagittal Saw	PR08200SB
Hall MicroFree Oscillating Saw	PR08300SB
Hall MicroFree Reciprocating Saw	PR08400SB
NEW Hall [®] MicroFree [™] Lithium Power Unit	
Hall MicroFree Lithium Power Unit	L3500SB
NEW Hall [®] MicroFree [™] Mini-Driver	
Hall MicroFree Mini-Driver	PR08500SB
NEW Hall [®] Small Lithium Battery	
Hall Small Lithium Battery	L3000SM
NEW Hall [®] MicroFree [™] Adapter and Charger	
Hall MicroFree Power Unit Adapter	L3500
Hall Lithium Battery Charger	L3000
NEW Hall [®] MicroFree [™] Sterilization Containers	
1/2 Size Small Bone Insert	PR08003T
1/2 Size Procedure Specific Insert	PR08001T
1/2 Size Rigid Container	TR12R
3/4 Size Rigid Container	TR34R
Accessories: Bur Guards	
Medium Bur Guard	1375-012
Long Bur Guard	1375-011
Extra-Long Bur Guard	1375-023

Local 727-392-6464CONMED.comToll Free 1-800-237-0169customer_service@CONMED.com