


Tonsillektomie, Laryngo-Bronchoskopie
Tonsillectomy, Laryngo-Bronchoscopy
Tonsilectomía, Laringo-Bronhoscopía
Tonsillectomie, Laryngo-Bronchoscopie
Tonsilectomia, Laringo-Bronhoscopia

39




Gaumenspalten-Instrumente

Cleft Palate Instruments


Instrumentos paladares

Instruments pour palais


Strumenti per palato




3,9 mm
39-001-18-07
Gaumenspalten-Häkchen, kurz
Cleft Palate Hooklet, short
Ganchito para el paladar, corto
Crochet pour palais, court
Uncino per palato, corto



5,5 mm
39-003-18-07
Gaumenspalten-Häkchen, lang
Cleft Palate Hooklet, long
Ganchito para el paladar, larga
Crochet pour palais, longue
Uncino per palato, lunga




3,3 mm
39-011-18-07
Gaumenspalten-Dissektor, kurz
Cleft Palate Dissector, short
Disecador para el paladar, corto
Dissecteur pour palais, court
Dissettore per palato, corto



3,7 mm
39-013-18-07
Gaumenspalten-Dissektor, lang
Cleft Palate Dissector, long
Disecador para el paladar, larga
Dissecteur pour palais, longue
Dissettore per palato, lunga




Salyer
39-001-18-07 -
39-027-18-07
18 cm / 7 1/8"



3,7 mm
39-021-18-07
Gaumenspalten-Elevatorium
Cleft Palate Instruments
Instrumentos paladares
Instruments pour palais
Strumenti per palato



3,6 mm
39-023-18-07
Gaumenspalten-Elevatorium
Cleft Palate Elevator
Elevador para el paladar
Élévateur pour palais
Elevatore per palato



2,9 mm
39-027-18-07
Gaumenspalten-Elevatorium, graduert
Cleft Palate Elevator, graduated
Elevador para el paladar, graduado
Élévateur pour palais, gradué
Elevatore per palato, graduato



4 mm
39-101-18-07
links
left
izquierdo
gauche
sinistra



4 mm
39-102-18-07
rechts
right
derecho
droit
destra

39-101-18-07 -
39-102-18-07
14 cm / 5 4/8"

Tonsillennmesser, Tonsillenelevatoren

Tonsil Knife, Tonsil Dissectors

Cuchillo, Disectores para amígdalas

Bistouri pour amygdales, Dissecteurs pour amygdales

Bisturi per tonsille, Disettori per tonsille



Tonsillektomie, Laryngo-Bronchoskopie

Tonsillectomy, Laryngo-Bronchoscopy

Tonsilectomía, Laringo-Bronchoscopia

Tonsillectomie, Laryngo-Bronchoscopie

Tonsilectomia, Laringo-Bronchoscopia



Fisher
39-220-22-07
22,5 cm / 8 7/8"

Brünings
39-228-23-07
23 cm / 9"



Hurd
39-281-22-07
22,5 cm / 8 7/8"
Tonsillenelevator / Retraktor
Tonsil Dissector / Retractor
Disector / Separador para amígdalas
Dissecteur / Rétracteur pour amygdales
Disettore / Retrattore per tonsille


Tonsillenelevatorien, Tamponstopfer

Tonsil Dissectors, Gauze Packers

Disectores para amígdalas, Taponeadores

Dissecteurs pour amygdales, Porte-tampon

Disettori per tonsille, Tamponatori





Tonsillenklemmen

Tonsil Forceps
Pinzas para amígdalas
Pincés pour amygdales
Pinze per tonsille

Tonsillektomie, Laryngo-Bronchoskopie

Tonsillectomy, Laryngo-Bronchoscopy
Tonsilectomía, Laringo-Bronchoscopia
Tonsillectomie, Laryngo-Bronchoscopie
Tonsilectomia, Laringo-Bronchoscopia



White
39-419-18-07
18 cm / 7 1/8"

Tonsillenfasszange
Tonsil Grasping Forceps
Pinza para amígdalas
Pince pour amygdales
Pinze per tonsille




Blohmke
39-423-20-07
20 cm / 7 7/8"

Tonsillenfasszange
Tonsil Grasping Forceps
Pinza para amígdalas
Pince pour amygdales
Pinze per tonsille



Schnidt
39-443-19-07 -
39-445-19-07



Schnidt
39-447-19-07 -
39-449-19-07



Tonsillenschnürer

Tonsil Snares

Cierra-nudos para amígdalas

Serre-noeuds pour amygdales

Serranodi per tonsille



Ø 0,3 mm
39-513-01-07

Schnürdraht, Packung mit 100 Stück
Snare Wire, package of 100 pieces
Alambre para asas, paquete con 100 piezas
Fil pour anses, paquet avec 100 pièces
Filo per legare, confezione di 100 pezzi

Tyding
39-525-01-07

mit einem geraden Aufsatz
with one straight tip
con una punta recta
avec une pointe droite
con una punta diritta



Broncho-Oesophagoskopiezangen

Broncho-Esophagoscopy Forceps

Pinzas para Broncho-Esofagoscopia

Pinces pour Broncho-Esofagoscopie

Pinze per Broncho-Esofagoscopia



Laryngo-Bronchoskopie

Tonsillectomy, Laryngo-Bronchoscopy


Tonsilectomía, Laringo-Bronchoscopía

Tonsillectomia, Laryngo-Bronchoscopie

Tonsilectomia, Laringo-Bronchoscopia




1/1



39-541-13-07

Schnitt
cut
corte
tranchant
tagliente
Ø 13 mm

Schaftlänge
shaft length
longitud del vástago
longueur de la tige
lunghezza dell'asta
120 mm



39-541-15-07

Schnitt
cut
corte
tranchant
tagliente
Ø 15 mm

Schaftlänge
shaft length
longitud del vástago
longueur de la tige
lunghezza dell'asta
120 mm

1/2

Hartmann
39-541-13-07 -
39-541-15-07

Tonsillenzange, durchschneidend
Tonsil Punchthrough-cutting
Pinzas para amígdalas, cortante
Pince pour amygdales, tranchante
Pinza per tonsille, tagliente


Broncho-Oesophagoskopiezangen

Broncho-Esophagoscopy Forceps
 Pinzas para Broncho-Esofagoscopia
 Pinces pour Broncho-Esofagoscopie
 Pinze per Broncho-Esofagoscopia



Brünings

39-584-20-07	22 cm / 8 5/8"
39-584-30-07	32 cm / 12 5/8"
39-584-40-07	41 cm / 16 1/8"



Brünings

39-587-20-07
 18 cm / 7 1/8"



Huber

39-581-01-07

für Verlängerungsrohre 39-584-20-07 - 39-587-20-07
 for interchangeable tubes 39-584-20-07 - 39-587-20-07
 para tubos de extension 39-584-20-07 - 39-587-20-07
 pour tubes rallonges 39-584-20-07 - 39-587-20-07
 per tubi allunghi 39-584-20-07 - 39-587-20-07

Auswechselbare Zangenaufsätze
Exchangeable Tips
Pinzas intercambiables
Embouts interchangeables
Pinze intercambiabili

Schnittbreite
Cutting width
Ancho de cortar
Largeur de coupe
Ampiezza di taglio




Fraenkel

39-600-01-07



Seiffert

39-604-01-07



Fraenkel

39-608-01-07



Brünings

39-616-02-07 - 39-616-03-07



Fig. 2
 39-616-02-07






Fig. 3
 39-616-03-07

Auswechselbare Zangenaufsätze

Exchangeable Tips
Pinzas intercambiables
Embouts interchangeables
Pinze intercambiabili

Laryngo-Bronchoskopie

Tonsillectomy, Laryngo-Bronchoscopy
Tonsilectomía, Laringo-Bronchoscopia
Tonsillectomie, Laryngo-Bronchoscopie
Tonsilectomia, Laringo-Bronchoscopia



Kehlkopfspiegel

Laryngeal Mirrors

Espejos laringeos

Miroirs laryngiens

Specchi laringei




Fig. 3	39-781-16-01	Ø 16 mm
Fig. 4	39-781-18-01	Ø 18 mm
Fig. 5	39-781-20-01	Ø 20 mm
Fig. 6	39-781-22-01	Ø 22 mm
Fig. 7	39-781-24-01	Ø 24 mm
Fig. 8	39-781-26-01	Ø 26 mm




Fig. 3	39-785-16-01	Ø 16 mm
Fig. 4	39-785-18-01	Ø 18 mm
Fig. 5	39-785-20-01	Ø 20 mm
Fig. 6	39-785-22-01	Ø 22 mm
Fig. 7	39-785-24-01	Ø 24 mm
Fig. 8	39-785-26-01	Ø 26 mm
Fig. 9	39-785-28-01	Ø 28 mm



1/2

rhodiumbeschichtet rhodium plated cubierto de rodio plaqué du rhodium coperto di rodio		
Fig. 00	39-782-08-04	Ø 8 mm
Fig. 0	39-782-10-04	Ø 10 mm
Fig. 1	39-782-12-04	Ø 12 mm
Fig. 2	39-782-14-04	Ø 14 mm
Fig. 3	39-782-16-04	Ø 16 mm
Fig. 4	39-782-18-04	Ø 18 mm
Fig. 5	39-782-20-04	Ø 20 mm
Fig. 6	39-782-22-04	Ø 22 mm
Fig. 7	39-782-24-04	Ø 24 mm
Fig. 8	39-782-26-04	Ø 26 mm
Fig. 9	39-782-28-04	Ø 28 mm
Fig. 10	39-782-30-04	Ø 30 mm

39-781-16-01 -
39-782-30-04

18,5 cm / 7 3/8"
plan, mit festem Handgriff
plane, with fixed handle
plano, con mango fijado
plat, avec manche fixe
plano, con manico fisso



1/2

rhodiumbeschichtet rhodium plated cubierto de rodio plaqué du rhodium coperto di rodio		
Fig. 00	39-786-08-04	Ø 8 mm
Fig. 0	39-786-10-04	Ø 10 mm
Fig. 1	39-786-12-04	Ø 12 mm
Fig. 2	39-786-14-04	Ø 14 mm
Fig. 3	39-786-16-04	Ø 16 mm
Fig. 4	39-786-18-04	Ø 18 mm
Fig. 5	39-786-20-04	Ø 20 mm
Fig. 6	39-786-22-04	Ø 22 mm
Fig. 7	39-786-24-04	Ø 24 mm
Fig. 8	39-786-26-04	Ø 26 mm
Fig. 9	39-786-28-04	Ø 28 mm
Fig. 10	39-786-30-04	Ø 30 mm

39-785-16-01 -
39-786-30-04


13,5 cm / 5 3/8"
plan, ohne Handgriff
plane, without handle
plano, sin mango
plat, sans manche
plano, senza manico

Kehlkopfspiegel -und Griffe

Laryngeal Mirrors, Laryngeal Mirror Handles
Espejos laringeos, Mangos para espejos laringeos
Miroirs laryngiens, Manche pour miroirs laryngiens
Specchi laringei, Manico per specchi laringei

Tonsillektomie, Laryngo-Bronchoskopie

Tonsillectomy, Laryngo-Bronchoscopy
Tonsilectomía, Laringo-Bronchoscopia
Tonsillectomie, Laryngo-Bronchoscopie
Tonsilectomia, Laringo-Bronchoscopia



Cr



1/1

Fig. 1
39-789-12-02
Ø 12 mm



Fig. 2
39-789-14-02
Ø 14 mm



Fig. 3
39-789-16-02
Ø 16 mm



Fig. 4
39-789-18-02
Ø 18 mm

1/2

Michel
39-789-12-02 -
39-789-18-02

27 cm / 10 5/8"

beweglicher Spiegel
moveable mirror
espejo movable
miroir mobile
specchio mobile



1/2

39-794-11-02
11 cm / 4 1/8"

➔ 39-785-08-01 -
39-786-30-04



1/2

39-798-12-02
12,5 cm / 4 7/8"

➔ 39-785-08-01 -
39-786-30-04

39

