

「Urofino」 Disposable Electronic Flexible Cystoscope



- Single use avoids cross infection and saves sterilization costs.
- Bigger working channel allows bigger instrument use.
- Taper-shape distal tip makes insertion easier.
- High resolution image makes procedure safer and easier.



Model	UV-100-E	UV-100T-E
Field of view	110°	110°
Depth of view	5-50mm	5-50mm
Deflection up/down	210°/130°	210°/130°
Insertion dia.	4.7mm	5.3mm
Distal end dia.	4.7mm	5.3mm
Working length	385mm	385mm
Working channel	2.2mm	2.6mm



Shanghai SeeGen Photoelectric Technology Co.,Ltd. Add.: 3 Floor, Building No.1 4299 Jindu Road, Minhang Dist., Shanghai China Email:endoexp@163.com

