

**NEW**

# Shortest Cycle Times in its Class – the New Belimed CLEANSTATION CS 750

Fast, economical and environmentally friendly – those are the characteristics of the new CLEANSTATION CS 750, the follow-up model of the Kleindienst/Belimed CLEANSTATION which was built over 1000 times. It is the right choice for efficient reprocessing of carts, containers, OR shoes, OR tables, beds and other items.



*CLEANSTATION CS 750 washer, disinfector and dryer for carts, containers, beds and other large volume washing items.*

## **Perfect quality with short cycle times**

With up to 10 cycles per hour, the CLEANSTATION CS 750 has an excellent throughput. The new, highly efficient jet spray system combined with sophisticated water piping and a flexible tank system produces excellent cleaning results within the shortest possible times, while energy and media consumption is minimized. High performance drying completely dries the washing items and allows fast reuse of equipment.

## **Easy and reliable handling**

The tinted glass pivot doors enable viewing of the illuminated washing chamber interior. The doors can be opened quickly and smoothly. Greater user comfort is provided by the operating panel which is uniform for all Belimed machines. It ensures clear intuitive menu navigation on an illuminated colour display. Important process data, such as remaining cycle time, readiness for loading or unloading as well as warnings are displayed by the patented process status display which is easily distinguishable from a distance.

**Belimed**  
Infection Control

System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed Austria: +43 3155 40699 0, Belgium: +32 10 42 02 40, China: +86 21 513 709 98, France: +33 3 89 63 65 40, Germany: +49 8631 9896 0, Hungary +36 1 318 86 97, Netherlands: +31 10 286 17 50, Slovenia: +386 1 7866 010, Switzerland: 0848 55 88 11, United Kingdom: +44 1403 738 811, USA: +1 843 216 7424, other countries: +41 41 449 78 88, [www.belimed.com](http://www.belimed.com)



Important process data such as the remaining cycle time, loading or removal readiness and error messages are displayed via our patented status display that is well visible from afar.

Illuminated Belimed touch panel CP-TOP with clear menu navigation.

The tilting floor of the CS 750 prevents water buildup on flat surfaces ensuring a fast drying time.

### Unique flexibility, optimal space utilisation

The available space is often limited, and this is where the CLEANSTATION CS 750 scores with an excellent ratio of usable chamber volume to footprint. Compact constructional design and a variety of available chamber lengths and widths enable implementation of a customized layout optimally suited to your requirements. To increase capacity, 2 to 5 additional modules can be integrated to form a multiple chamber system with conveyor system.

### Easy-to-service

Right from the start, the CLEANSTATION CS 750 was designed for easy servicing and maintenance. All units and the control system are clearly arranged and housed in a separate servicing space optionally located on the right or left for easy accessibility.

## Technical data

### Dimensions

Type	CS 750S	CS 750M	CS 750L
Usable chamber dimensions H × W × D (mm)	2000 × 1000 × 1350	2000 × 1000 × 2300	2000 × 1000 × 2950
External dimensions H × W × D (mm)	3010 × 2500 × 2250	3010 × 2500 × 3200	3010 × 2500 × 3850
Loading height	floor level	floor level	floor level
Pit depth (mm)	120	120	120

Capacity per cycle	Pieces
Carts	1–3
Containers	12
OR shoes	80 (40 pairs)
OR tables	1
OR dishes ø 200–300 mm	24
Beds	1

**Belimed**  
Infection Control

### System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed Austria: +43 3155 40699 0, Belgium: +32 10 42 02 40, China: +86 21 513 709 98, France: +33 3 89 63 65 40, Germany: +49 8631 9896 0, Hungary +36 1 318 86 97, Netherlands: +31 10 286 17 50, Slovenia: +386 1 7866 010, Switzerland: 0848 55 88 11, United Kingdom: +44 1403 738 811, USA: +1 843 216 7424, other countries: +41 41 449 78 88, [www.belimed.com](http://www.belimed.com)